



## Experience, Quality and Technology

With 45+ years experience, we are constantly working toward better solutions...

**medbar<sup>®</sup>**

 +90 530 938 8175



STERILE EO CE 2292  

## Arterial Cannula

Ref: 304

*304 01 Arterial Cannula 20G    304 02 Arterial Cannula 18G    304 03 Arterial Cannula 22G*

- Atraumatic needle tip: specifically designed to minimize injurious effects or trauma.
- 18 G – 1.30 mm x 45 mm - 80 ml/min.    20 G – 1.10 mm x 45 mm - 50 ml/min.    22 G – 0.90 mm x 33 mm - 30 ml/min.
- Easily fixed into place with soft wings.
- Built in air embolism protection: Blood back-flow and blood contamination are also avoided with the on/off switch.
- The flow-block mechanism can be used easily with one hand.

**Alternative Names** Arterial Catheter • Radial Cannula • Femoral Cannula / Catheter

## INFUSION



STERILE EO CE 2292  

## IV Flow Controller

Ref: 233 01

- Gravity type flow controller.
- Scale range from 2 to 350 ml/hour.
- Constant flow rate.
- Universal male and female luer connections.
- Y-port injection.
- Standard length tubing 45cm.
- Different length tubing available.

**Alternative Names** IV Flow Regulator • IV Precision Flow Device



STERILE EO CE 2292

## Arthroscopy Set

Ref: 203

*203 01 Arthroscopy Set with Manual Pressure Pump*   *203 02 Arthroscopy Set*

- Y-type irrigation set.
- Non-vented spike.
- Airproof incision ends.
- Clamps.
- Y-connector.
- Air-freeing lid.

**Alternative Names** Arthroscopy Fluid Management Device (Irrigation Console Accessory)

## SURGICAL PRODUCTS



STERILE EO CE 2292

## Skin Marking Set

Ref: 241

*241 01 Skin Marking Set (Fine Tip)*   *241 03 Skin Marking Set (Ultra Fine Tip)*

- Fine tip is considered the traditional surgical tip.
- Ultra fine tip is widely used for plastic and facial surgery applications.
- Biocompatible ink.
- Flexible ruler comes in cms and inches.
- Insoluble.
- Non-toxic and non-tattooing.

**Alternative Names** Surgical Skin Marking Set • Surgical Marker



BIOCIDAL  

## Impregnated Surgical Hand Brush

Ref: 205

*205 01 Single Use Surgical Hand Brush Impregnated with 4% Chlorhexidine Gluconate*

*205 02 Single Use Surgical Hand Brush Impregnated with 7.5% Providone-Iodine*

- Ergonomic brush design (easy to use).
- Easy open package.
- Soft bristles to avoid skin irritation.
- With or without nail pick.
- Smooth and soft sponge absorbs up to 90 ml of liquid.
- Chlorhexidine Gluconate and Providone-Iodine options.
- For external use only.

### Alternative Names

Impregnated Scrub Brush • Surgical Brush • Impregnated Surgical Brush • Surgical Sponge



STERILE EO  

## Dry Sterile Surgical Hand Brush

Ref: 205 03

- Ergonomic brush design (48 x 79 x 34mm).
- Easy open package.
- Soft bristles for no skin irritation.
- With or without nail pick.
- Smooth and soft sponge absorbs up to 90 ml liquid.
- For external use only.

### Alternative Names

Disposable Surgical Brush • Sterile Surgical Brush • Sterile Dry Surgical Brush

## Karman Cannula Manual Vacuum Aspirator Double Valve

Ref: 209

*209 02 Karman Cannula Manual Vacuum Aspirator Double Valve Procedure Pack (Aspirator + Cannulas)*

*209 05 Karman Cannula Manual Vacuum Aspirator Double Valve (No Cannulas)*



STERILE EO CE 2292

- Includes extra O-ring, extra collar stop and 2cc silicone oil.
- Adaptors included.
- Stable vacuum: each product's vacuum is individually tested.
- Portable.
- 60cc barrel capable of producing a vacuum of 609.6-660.4 mm/Hg.
- Safe and effective.
- Low cost.
- Appropriate for many different clinical settings.
- High patient and provider satisfaction.
- Outpatient Procedure.
- Products of Conception (POC) easily visible.
- Single use device.
- Our Karman Cannula Manual Vacuum Aspirator Double Valve Procedure Pack includes sterile Karman Cannulas sizes 6, 7, 8, 9, 10 and 12.

### Alternative Names

Suction Aspirator • Miscarriage Management Kit • Manual Uterine Aspirator • Uterine Evacuation Service • Gynecological Aspirator • Endometrial Aspirator • Aspirating Syringe • Aspirator Kit



STERILE EO CE 2292

### Karman Cannula Manual Vacuum Aspirator Single Valve

Ref: 209

**209 01** Karman Cannula Manual Vacuum Aspirator Single Valve Procedure Pack (Aspirator + Cannulas)

**209 04** Karman Cannula Manual Vacuum Aspirator Single Valve (No Cannulas)

- Includes extra O-ring, extra collar stop and 2cc silicone oil.
- Stable vacuum: Each product's vacuum is individually tested.
- Portable.
- 60cc barrel capable of producing a vacuum of 609.6-660.4 mm/Hg.
- Safe and effective.
- Low cost.
- Appropriate for many different clinical settings.
- High patient and provider satisfaction.
- Outpatient Procedure.
- Single use device.
- Our Karman Cannula Manual Vacuum Aspirator Single Valve Procedure Pack includes sterile Karman Cannulas sizes 4, 5 and 6.

#### Alternative Names

Suction Aspirator • Miscarriage Management Kit • Manual Uterine Aspirator • Uterine Evacuation Service • Gynecological Aspirator • Endometrial Aspirator • Aspirating Syringe • Aspirator Kit

### Karman Cannula

Ref: 208

**208 03** No: 3 **208 04** No: 4 **208 05** No: 5 **208 06** No: 6 **208 07** No: 7

**208 08** No: 8 **208 09** No: 9 **208 10** No: 10 **208 12** No: 12

- Total Length: 22 cm.
- Smooth catheter surface.
- Rounded tip.
- Individually packed in peelable pouch pack.
- 3-12 mm diameter options.
- Matching color codes for adaptor for easy pairing (when necessary).

#### Alternative Names

MVA Cannula • Intra-uterine Suction Cannula



STERILE EO CE 2292

## Endometrial Suction Curette

Ref: 232 01



STERILE EO CE 2292

- Outside diameter: 3 mm.
- Length (with pistons): 26 cm.
- Appropriate for endometrial biopsy and for use in menstrual extractions.
- Smooth catheter surface.
- Flexible and soft, allowing the curette to adapt to each unique uterine cavity.
- Has four holes on tip.
- Dilation not required due to the small size of the tip.
- Provides for a simple and atraumatic procedure.
- Procedure carried out in less than one minute.

### Alternative Names

Biopsy Cannula • Endometrial Biopsy Cannula • Intra-uterine Suction Cannula

## Endometrial Suction Curette with Syringe

Ref: 232 02



STERILE EO CE 2292

- Outside diameter: 3mm.
- Appropriate for endometrial biopsy and for use in menstrual extraction.
- Smooth catheter surface.
- Provides for a simple and atraumatic procedure.
- Easier sampling with the syringe.

### Alternative Names

Endometrial Aspiration Cannula with Syringe • Intra-uterine Aspiration Cannula with Syringe • Endometrial Biopsy Cannula with Syringe



### Insemination Cannula

Ref: 286 01

- The sliding handle on the cannula helps to adjust to the true distance of the uterus.
- The even distribution of sperm is ensured by the openings on both sides of the tip.
- Atraumatic uterine and cervical penetration.
- Universal connector to connect the tube to a syringe.
- The luer lock on the connector allows for secure locking to take place.

#### Alternative Names

Intra-uterine Insemination Cannula • Embryo Transfer Cannula



### Insemination Cannula with Guided Wire

Ref: 286 02

- With the aid of the guided wire the catheter can form to the unique shape of each individual uterus.
- The sliding handle on the cannula helps to adjust to the true depth of the uterus.
- Atraumatic uterine and cervical penetration.
- Universal connector to connect the tube to a syringe.
- The luer lock on the connector allows for secure locking to take place.

#### Alternative Names

Intra-uterine Insemination Cannula with Guided Wire • Embryo Transfer Cannula with Guided Wire





STERILE EO CE 2292  

## Cervix Brush

Ref: 211

*211 01 Cervix Brush (Sterile) 211 02 Cervix Brush (Non-Sterile) 211 03 Cervical Brush (Sterile)*

- Designed to collect exfoliated cervical cells.
- The harvested cells should be processed with normal or liquid based cytology (LBC).
- In order to collect as many cells as possible the brush has been designed with soft, flexible, smooth and dense bristles.
- Total length: 20 cm.
- If desired, the brush can be easily separated from the brush handle.
- Sterile and non-sterile options.

**Alternative Names** Cytology Brush • Endocervical Brush



STERILE EO CE 2292  

## Pouch Perforator

Ref: 210 01

- The specially designed tip allows the amniotic membrane (pouch) to be perforated in a single motion.
- Rounded tip designed to rupture the amniotic membrane without injury to the mother or fetus.
- Total length 26 cm.

**Alternative Names** Amnihook • Amniotic Pouch Perforator • Amnio Hook • Amniotic Membrane Perforator/Breaker



STERILE EO CE 2292  

## Smear Spatula

Ref: 240 01

- Provides easy and smooth excavation from the cervix mouth without causing harm to the patient.
- One side is for the purpose of excavation and the other side is for the purpose of transferring.
- Extended tip better for collecting endocervical cells.

**Alternative Names** Pap Smear Spatula • Cervical Spatula • Ayre's Spatula • Aylesbury Spatula



### Cardboard Camera Cover

Ref: 202 01

- Centralized hanger hole.
- Telescopically folded.
- Single layer thickness is 40 micron.
- Eight layer folding allows the cover to pass over the camera easier and faster.
- Compatible with all camera models.
- Different sizes are available upon request.

**Alternative Names** Disposable Endoscopy Camera Cover • Sterile Endoscopy Camera Cover



### Video Camera Cover with Elastic Tip

Ref: 202 07

- Specially designed elastic tip provides easy application.
- Telescopically folded.
- Single layer thickness is 40 micron.
- Eight layer folding allows the cover to pass over the camera easier and faster.
- Compatible with all camera models.
- Different sizes are available upon request.

**Alternative Names** Disposable Endoscopy Camera Cover • Sterile Endoscopy Camera Cover



### Circled Ring Camera Cover

Ref: 202 02

- 14 x 250 cm.
- Perfectly folded for fast application.
- Has bigger mid-hole than any other Camera Cover on the market.
- Easy to fix adhesive straps (extra strap provided with each product).
- Different sizes are available upon request.

**Alternative Names** Disposable Endoscopy Ring Camera Cover



STERILE EO CE 2292

### Circled Ring Camera Cover with Elastic Tip

Ref: 202 09

- Specially designed elastic tip provides easy application.
- 14 x 250 cm.
- Perfectly folded for fast application.
- Has bigger mid-hole than any other Camera Cover on the market.
- Easy to fix adhesive strap.
- Different sizes are available upon request.



STERILE EO CE 2292

### Microscope Drape

Ref: 204

*204 01 Microscope Drape (48 mm) 204 02 Microscope Drape (65 mm) 204 03 Microscope Drape (68 mm)  
204 07 Microscope Drape (65 mm Sloped Lens)*

- Specially folded for easy application.
- Super transparent lens system.
- Adhesive straps secure the drape safely in place.
- Accessories are easy to use when the drape is in place.
- Compatible with Carl Zeiss, Leica and Moller-Wedel operating microscope models.
- Different sizes are available upon request.



STERILE EO CE 2292

### C-Arm Drape

Ref: 227

*227 04 C-Arm Drapes with Clips 227 08 C-Arm Drape with Clips (Medium) 211 11 C-Arm Drape with Clips (Small)*

- Includes two Medical Camera Drapes and one C-Arm Drape with clips.
- The camera arm is easily covered by the clip-on cover and prevents contamination.
- Locking clips are included.
- No contamination: perfectly surrounds the equipment.
- Compatible with Siemens, Philips, Ziehm and GE C-Arm machines.



STERILE EO CE 2292

## Fluoroscopy Drape

Ref: 222 01

- 90 x 60 cm standard size.
- Compatible with all devices.
- Different sizes are available upon request.

Alternative Names Banded Bag



STERILE EO CE 2292

## Scopy Drape

Ref: 227 01

- Standard 60 cm diameter.
- Transparent: does not inhibit the functioning of the device during use.
- Different sizes are available upon request.

Alternative Names Dome Bag

# ENDOSCOPY



STERILE EO CE 2292

## Endoscopy Mouthpiece

Ref: 207

- 207 01 Endoscopy Mouthpiece with Strap (Adult)*
- 207 02 Endoscopy Mouthpiece with Strap (Pediatric)*
- 207 03 Endoscopy Mouthpiece without Strap*
- 207 04 Endoscopy Mouthpiece with Oxygen Entry Port*
- 207 06 Endoscopy Mouthpiece without Strap and Tooth Protector*
- 207 07 Endoscopy Mouthpiece with Oxygen Entry and Strap*

- Smooth edges and smooth passing surface.
- Silicone tooth protector.
- Strap to secure comfortably in place during the procedure.
- Compatible with mouth anatomy.
- Fits different sizes of endoscopes.
- Endoscopy Mouthpiece with Oxygen Entry Port allows oxygen to be given through the oxygen entrance.

Alternative Names Endoscopy Mouth Guard • Endoscopy Bite Block



STERILE EO CE 2292

## Endo Seal Bag

Ref: 279

279 01 Endoseal Bag (Medium)    279 02 Endoseal Bag (Large)

- Patent protected.
- Special sterile endoscope storage and reversible carrying bag (and liner).
- Transparent window to indentify contents easily.
- Clean/contaminated markings to easily distinguish as clean or dirty.
- Completely contains the soiled equipment, reducing the risk of spreading infection.
- Sealing of the bag prevents contamination.



CE

## Endo Clean Sponge

Ref: 290 01

- Used for pre-cleaning of endoscopes to remove gross contaminants.
- Made of special elastic foam for a soft touch when cleaning sensitive endoscopes.
- The mid hole fits the outer sheath of endoscopes perfectly.
- Safe for use on all scopes (flexible/rigid) and instruments.



STERILE EO CE 2292

## Phototherapy Eye Band

Ref: 201

201 01 Phototherapy Eye Band (Small)    201 02 Phototherapy Eye Band (Medium)  
201 03 Phototherapy Eye Band (Large)

- Special anti-perspirant fabric.
- Blocks 100% of light.
- Adjustable for a perfect fit.

**Alternative Names** Phototherapy Ocular Protector • Phototherapy Mask



### Phototherapy Eye Band (Y-band)

Ref: 236

236 01 Y-band (Small) 236 02 Y-band (Medium) 236 03 Y-band (Large)

- Unique Y-shaped design.
- Blocks 100% of light.
- Elastic for a perfect fit.
- Eye-pad material blocks light from reaching the eyes, while headband allows phototherapy light to reach the baby's head.
- Provides eye protection for neonatal phototherapy.



### Valve Emesis Bag

Ref: 228 01

- 1000 ml capacity.
- Graduated scale allows for the precise monitoring of fluid loss.
- Advanced valve system (no spills, no smell, and no leakage).
- Advanced anti-reflex design.

**Alternative Names** Valve Vomit Bag



### Emesis Bag

Ref: 228 02

- Includes an ultra absorbent pad that absorbs up to 450 ml of liquid.
- Constriction system with double-layered rope that allows the bag to be securely closed.

**Alternative Names** Vomit Bag



STERILE | EO | CE 2292 | MDP | X

## Limb Holder

Ref: 238

*238 01 Limb Holder (Adult/Hand) 238 02 Limb Holder (Adult/Ankle)*  
*238 03 Limb Holder (Infant) 238 04 Limb Holder (Children)*

- Perfect for emergency rooms.
- Secure tying system.
- Does not irritate skin.
- Individually packaged in peelable pouch pack.

**Alternative Names** Limb Restraint



STERILE | EO | CE 2292 | MDP | X

## Tracheostomy Tube Fixer

Ref: 262

*262 01 Tracheostomy Tube Fixer (Small) 262 02 Tracheostomy Tube Fixer (Medium)*  
*262 03 Tracheostomy Tube Fixer (Large)*

- 34 x 7 x 1 cm.
- Special fabric allows oxygen and steam to pass through.
- The surface that comes into contact with the skin is made of special anti-perspirant fabric.
- Easily adjustable to the desired length.

**Alternative Names** Tracheostomy Tube Holder • Tracheostomy Tube Fastner



STERILE | EO | CE 2292 | MDP | X

## Endotracheal Tube Fixer

Ref: 281

*281 01 Endotracheal Tube Fixer (Small) 281 02 Endotracheal Tube Fixer (Medium)*  
*281 03 Endotracheal Tube Fixer (Large)*

- Special fabric allows oxygen and steam to pass through.
- The surface that comes into contact with the skin is made of special anti-perspirant fabric.
- Easily adjustable to the desired length.

**Alternative Names** Endotracheal Tube Holder • Endotracheal Tube Fastner



### Spirometer Filter

Ref: 273

273 01 SFT 33 mm 273 02 SFT 30mm 273 03 SFT 33 mm with Latch 273 04 SFT 30 mm with Latch 273 05 SFT 26 mm 273 06 SFT 26 mm with Latch  
273 07 SFT 33 mm with Latch (Large) 273 08 SFT 30 mm with Latch (Large) 273 09 SFT 33 mm (Large) 273 10 SFT 30 mm (Large)

- Made out of electrostatic material.
- Bacterial Filter Efficiency Test (BFE) 99.99%.
- Viral Filtration Efficiency Test (VFE) 99.99%.
- Easy to open OPP packaging.
- Dead space less than 80 ml.
- Manufactured without the use of chemical adhesives.

**Alternative Names** Spirometry Mouthpiece



### Mucous Aspirator

Ref: 229

229 01 Mucous Aspirator 40 ml 229 02 Mucous Aspirator 20 ml 229 03 Mucous Aspirator 100 ml 229 04 Mucous Aspirator 10 ml

- Compatible with all suction catheters between 4-20 ch.
- A sticker is included in the set to use for patient information.
- Extra covers are available.

**Alternative Names** Mucous Sucker • Mucous Extractor • Mucous Extraction Container  
Mucous Trap • Mucous Suction Pump



### Mucous Aspirator with Silicone Tube

Ref: 305 01

- After the procedure the container can be easily sealed via the use of the silicon tube.
- Compatible with all suction catheters between 4-20 ch.
- Up to 40 ml. of aspiration capacity.
- Clearly defined accurate scale.

**Alternative Names** Mucous Extractor with Silicone Tube • Mucous Trap with Silicone Tube • Mucus Specimen Trap  
Suction Trap • Lukens Trap





CE 2292

## Valve Urine Bag

Ref: 225 03

*225 01 Valve Urine Bag*      *225 03 Valve Urine Bag with Discharger*

- 1500 ml Scale: The amount of urine can be easily monitored.
- Advanced valve system (no spills, no leakage).
- Discharge system for easy emptying (Optional).
- Anatomically designed for use by both males and females.
- The special handle system allows for the bag to be hung anywhere easily.



CE



## Easy Flush

Ref: 239 01

- Ultra absorbent pad.
- Up to 450 ml liquid absorption capacity.
- It includes a rope that allows the bag to be closed.

Alternative Names      Bedpan Liner

# medbar®

Medbar Tıbbi Malzemeler Turizm San. Tic. A.Ş.  
1142 Sk. No:35 Sarnıç Gaziemir - İzmir / Türkiye  
T.: +90 232 281 60 03 F.: +90 232 281 66 48

[www.medbar.com.tr](http://www.medbar.com.tr) [info@medbar.com.tr](mailto:info@medbar.com.tr)

+90 530 938 81 75

