



Endometrial Suction Curette with Injector



REF 232 02



STERILE EO

Endometrial Suction Curette with Injector Features

- * Outside diameter: 3 mm.
- * Appropriate for endometrial biopsy and for use in menstrual extractions.
- * Smooth catheter surface.
- * Easy to use with the special design.
- * Allows a simple and non-traumatized application.
- * Easier sampling with the injector.

Quality

- * Manufactured under ISO 13485:2012 quality management standard.
- * CE certified.

Bio-Compatibility

- * Ethylene-oxide sterilized.
- * Latex free.

Shelf Life

- * 5 years.



REF NO	PRODUCT NAME	FEATURES	STERILIZATION	PACKAGING UNITS
23202	Endometrial Suction Curette with Injector	<ul style="list-style-type: none"> • Outside diameter: 3 mm. • Appropriate for endometrial biopsy and for use in menstrual extractions. • Smooth catheter surface. • Easy to use with the special design. • Allows a simple and non-traumatized application. • Easier sampling with the injector. 	ETO	50 pcs / inner bag 300 pcs / outer box Box Dimensions (cm): 42*58*37 Box Weight: 7,00 kg