

medbar[®]

C-Arm Drape with Clips



REF 227 04/08/11



CE 2292

STERILE EO



C-Arm Drape with Clips

Intended Use

- C-Arm Cover is used for isolating the contamination source to prevent contamination risk coming from hospital devices to media, patient or personnel during the surgical operations from contamination.
- C-Arm Covers are sterile and disposable medical devices with transient use.
- The devices are intended for single use.
- Main performance indicators for C-Arm Covers are sterility, to cover the device without preventing device use and to provide sterile media by preventing the contamination from the device to media.

Components of The Product

Medbar brand C-Arm Drape with Clips consists of :

- 2 piece scopy/fluoroscopy drapes,
- One Drape with 6 clips.

Product Features

- Specially folded for easy application.
- Easy to use directional sticker labels for correct application.
- Two durable straps to properly secure the device to the machine.
- More than different sizes available.
- Customized sizes are available upon request.



Quality

- CE marked according to 93/42/EEC Medical Device Directive and classified as Class I sterile medical device.
- Manufactured under ISO 13485:2012 Quality Management System standard.
- Notified Body: UDEM LTD A.Ş. (2292)
- EC Certification No: M.2016.106.7000
- GMDN Code: 43970

Sterilization

- Sterilized with Ethylene Oxide.
- Validated for ISO 11135 standard.

Shelf Life

- 5 years.

REF NO	PRODUCT NAME	STERILIZATION	CLIPS
227-04	C-Arm Drape 75*90 cm (Big)	EO	16*13*13 cm
227-08	C-Arm Drape 75*75 cm (Medium)	EO	13*10*10 cm
227-11	C-Arm Drape 75*75 cm (Small)	EO	13*7*10 cm

C-Arm Cover with Clips

Technical Drawing

C-Arm Cover with Clips

