

medbar[®]

Dry Sterile Surgical Hand Brush



REF 205 03



STERILE EO

BIOCIDAL

Dry Sterile Surgical Hand Brush

Intended Use

- Used to clean hands before surgeries.

Components of The Product

Medbar brand Dry Hand Brush consists of :

- Brush,
- Sponge,
- Nail Pick.

Product Features

- Ergonomic brush design (48*79*34 mm).
- Easy open package.
- Soft bristles to avoid skin irritation.
- With or Without nail prick.
- Smooth and soft sponge absorbs up to 90 ml liquid.
- Chlorhexidine Gloconate and Povidine-Iodine options.
- For external use only.

Quality

- Sterile handbrush is not classified as medical device according to European Commission.
- Manufactured under ISO 13485:2012 Quality Management System standard.

Sterilization


- Sterilized with Ethylene Oxide.

Shelf Life

- 5 years.

Biocompatibility

- Latex-free.



REF NO	PRODUCT NAME	STERILIZATION	PACKAGING
			40 pcs inner box/400 pcs outer box
205 03	Dry Sterile Surgical Hand Brush	EO	Box Dimensions: 67*35*54 cm
			Box Weight: 8,80Kg



Medbar Tıbbi Malzemeler Turizm San. Tic. A.Ş.
1142 Sk. No:35 Sarnıç Gaziemir-İzmir/Turkey
T: +90 232 2816003 - F: +90 232 2816648
www.medbar.com.tr - info@medbar.com.tr