

# medbar<sup>®</sup>

## Impregnated Hand Brushes



REF 205



STERILE EO

BIOCIDAL

## Impregnated Hand Brushes

### Intended Use

- Used to clean hands before surgeries.

### Components of The Product

Medbar brand Impregnated Hand Brushes consist of :

- Brush,
- Sponge,
- Biocidal.

### Product Features

- Ergonomic brush design (48\*79\*34 mm).
- Easy open package.
- Soft bristles to avoid skin irritation.
- With or Without nail prick.
- Smooth and soft sponge absorbs up to 90 ml liquid.
- Chlorhexidine Gloconate and Povidine-Iodine options.
- For external use only.

### Quality

- Biocidal Licensed.(Impregnated).
- Manufactured under ISO 13485:2012 Quality Management System standard.

### Sterilization

- Non-sterile.

### Shelf Life

- 2 years.

### Biocompatibility

- Latex-free.



REF NO	PRODUCT NAME	STERILIZATION	PACKAGING
205 01	Single Use Surgical Hand Brush 4% Chlorhexidine	-	40 pcs inner box/400 pcs outer box
205 02	Single Use Surgical Hand Brush 7.5% Povidine-Iodine	-	Box Dimensions: 67*35*54 cm



Medbar Tıbbi Malzemeler Turizm San. Tic. A.Ş.  
1142 Sk. No:35 Sarnıç Gaziemir-İzmir/Turkey  
T: +90 232 2816003 - F: +90 232 2816648  
www.medbar.com.tr - info@medbar.com.tr