

# medbar®

## Smear Spatula



REF 240 01



CE 2292

STERILE EO



## Smear Spatula

### Intended Use

- The device is used in the collection of a pap smear screening procedure for cervical cancer.
- Smear Spatula are used by professional users generally for taking sample from the vagina for smear test.
- Smear Spatula is a sterile, single-use medical device with transient use.
- Main performance indicators for Smear Spatula are sterility, ability to collect the sample, allow to transmit sample to microscope slide and not to give any discomfort to patient during the sampling.

### Raw Material

- Polyoxymethylene.

### Quality

- Manufactured under ISO 13485:2012 Quality Management System standard.
- CE marked according to 93/42/EEC Medical Device Directive and classified as Class I sterile medical device.
- Notified Body: UDEM Ltd. Şti. (2292)
- EC Certification No: M.2016.106.7000
- GMDN Code: 32679
- Barcode: 8698881910758
- Ref No: 240 01

### Biocompatibility


- Latex free.

### Sterilization

- Sterilized with Ethylene Oxide.
- Validated with ISO 11135 standard.

### Shelf Life

- 5 years.



REF NO	PRODUCT NAME	STERILIZATION
240 01	Smear Spatula	EO



Medbar Tıbbi Malzemeler Turizm San. Tic. A.Ş.  
1142 Sk. No:35 Sarnıç Gaziemir-İzmir/Turkey  
T: +90 232 2816003 - F: +90 232 2816648  
www.medbar.com.tr - info@medbar.com.tr