

medbar[®]

Endo Clean Sponge



REF 290 01



Endo Clean Sponge

Product Features

- Used for pre-cleaning of endoscopes to remove gross contaminants.
- Made of special elastic foam for a soft touch when cleaning sensitive endoscopes.
- The mid hole fits the outer sheath of endoscopes perfectly.
- Safe for use on all scopes (flexible/rigid) and instruments.

Quality

- CE marked according to 93/42/EEC Medical Device Directive and classified as Class I sterile medical device.
- Manufactured under ISO 13485:2012 Quality Management System standard.

Biocompatibility

- Biologically evaluated according to ISO 10993-1.
- The product does not incorporate a substance of human blood derivative, animal originated tissues, phthalates, medicinal product, latex, radioactive material and electromagnetic waves.

Sterilization

- Endoclean Sponge is a non-sterile product.

Shelf Life

- 5 years.

