

# medbar<sup>®</sup>

## Valve Urine Bag with Discharger



REF 225 03



CE 2292



## Valve Urine Bag with Discharger

### Intended Use

- Valve Urine Bag is used for collecting urine in a hygienic media and measure the urine volume when necessary.
- Urine is collected into a scaled plastic bag, and unidirectional valve secures the urine to remain within the Valve Urine Bag and prevents the urine or odour to leak out.

### Components of The Product

Medbar brand Valve Urine Bag with Discharger consists of :

- Bag,
- Discharger sealed to the bag,
- Discharger stopper,
- Circle sealed to bag,
- Handle parts.

### Product Features

- 1500 ml scale: the amount of the urine can be easily monitored.
- Advanced valve system ( no spills, no leakage).
- Discharge system for easy emptying (Optional).
- Anatomically designed for use by both males and females.
- The special handle system allows for the bag to be hung anywhere easily.
- The discharger on the Valve Urine Bag with Discharger is intended for discharging the urine out.

### Quality

- CE Certified.
- Manufactured under ISO 13485:2012 Quality Management System standard.

### Sterilization

- Non-sterile.

### Shelf Life



- 5 years.

### Biocompatibility

- Phthalate free.
- Latex free.

### Warnings

- Re-use of the products may lead to infection or other illnesses .



REF NO	PRODUCT NAME	STERILIZATION
225 01	Valve Urine Bag	-
225 03	Valve Urine Bag with Discharger	-

