

# medbar®

**Medbar Tıbbi Malzemeler Turizm San. Tic. A.Ş.**  
1142 Sk. No:35 Sarnıç Gazıemir - İzmir / Türkiye

T: +90 232 281 60 03  
F: +90 232 281 66 48  
info@medbar.com.tr

[www.medbar.com.tr](http://www.medbar.com.tr)



**Experience, Quality and Technology**

With 45 + years experience, we are constantly working toward better solutions.

# medbar®

info@medbar.com.tr  
www.medbar.com.tr

As of today, Medbar has more than 200 employees and exports to more than 70 countries in the world.

We aim at major projects with our vision of  
**“Enhancing the quality of life by maximizing our contribution to the healthcare process”.**





STERILE EO CE 1984

**CYLINDRICAL IV FLOW CONTROLLER**

Ref: 133

- 133 01 Cylindrical IV Flow Controller
- 133 02 Cylindrical IV Flow Controller Long (62cm)
- 133 03 Cylindrical IV Flow Controller Rotary Luer Lock
- 133 04 Cylindrical IV Flow Controller without Y-Port
- 133 05 Cylindrical IV Flow Controller Needle Free

- It works on the principle of gravity.
- 5 to 250ml/h scale range is available. (Different ranges between 5 and 250ml/h can be provided upon request.)
- It has accurate and constant flow rate.
- Universal male and female luer connections are available.
- Needle and needle-free y-port injection options are available.
- Standard set length is 47cm. (Different hose lengths are available.)

**Alternative Names**

- IV Flow Controller
- IV Precision Flow Device
- Controflow

**IV FLOW CONTROLLER**


Ref: 233

- 233 01 IV Flow Controller with Y-Port
- 233 04 IV Flow Controller without Y-Port
- 233 05 IV Flow Controller Needle Free

- It works on the principle of gravity.
- 2 to 350 ml/hour scale range is available.
- It has accurate and constant flow rate.
- Universal male and female luer connections are available.
- Needle and needle-free y-port injection options are available.
- Standard hose length is 47 cm.(Different hose lengths are available.)

**Alternative Names**

- IV Flow Regulator
- IV Precision Flow Device
- Controflow



STERILE EO CE 1984



STERILE EO CE 1984

**EXTENSION LINE (250 PSI)**

Ref: 294

- 294 01 Extension Line 30cm (250 PSI)
- 294 02 Extension Line 50cm (250 PSI)
- 294 03 Extension Line 60cm (250 PSI)
- 294 04 Extension Line 75 cm (250 PSI)
- 294 05 Extension Line 90 cm (250 PSI)
- 294 06 Extension Line 100 cm (250 PSI)
- 294 07 Extension Line 120cm (250 PSI)
- 294 08 Extension Line 150cm (250 PSI)

- 2.70 mm-4.10 mm standard diameter.
- Universal male and female luer connectors are available.
- There are varieties with M/F and M/M connectors in sizes 30cm, 50cm, 60cm, 75cm, 90cm, 100cm, 120cm, 150cm.
- It is resistant to pressure up to 250 PSI (17 bar).



STERILE EO CE 1984

**PRESSURE RESISTANT EXTENSION LINE (800 PSI)**

Ref: 295

- 295 01 Pressure Resistant Extension Line 30 cm (800 PSI)
- 295 02 Pressure Resistant Extension Line 50 cm (800 PSI)
- 295 03 Pressure Resistant Extension Line 60 cm (800 PSI)
- 295 04 Pressure Resistant Extension Line 75 cm (800 PSI)
- 295 05 Pressure Resistant Extension Line 90 cm (800 PSI)
- 295 06 Pressure Resistant Extension Line 100 cm (800 PSI)
- 295 07 Pressure Resistant Extension Line 120 cm (800 PSI)
- 295 08 Pressure Resistant Extension Line 150 cm (800 PSI)

- 1.50 mm-3.00 mm standard diameter.
- Appropriate for all types of infusion treatment.
- Universal male and female luer connectors are available.
- There are varieties with M/F and M/M connectors in sizes 30cm, 50cm, 60cm, 75cm, 90cm, 100cm, 120cm, 150cm.
- It is resistant to pressure up to 800 PSI (55 bar).



STERILE EO CE 1984

**PRESSURE RESISTANT EXTENSION LINE (1200 PSI)**

Ref: 297

- 297 01 Pressure Resistant Extension Line 100 cm (1200 PSI)
- 297 02 Pressure Resistant Extension Line 120 cm (1200 PSI)
- 297 03 Pressure Resistant Extension Line 150 cm (1200 PSI)
- 297 04 Pressure Resistant Extension Line 180 cm (1200 PSI)
- 297 05 Pressure Resistant Extension Line 200 cm (1200 PSI)
- 297 06 Pressure Resistant Extension Line 250 cm (1200 PSI)
- 297 07 Pressure Resistant Extension Line 300 cm (1200 PSI)

- 1.50 mm-4.10 mm standard diameter.
- Various length options are available between 100-300 cm.
- Appropriate for all types of infusion treatment.
- Universal male and female luer connectors are available.
- It is resistant to pressure up to 1200 PSI (83 bar).

**SINGLE VALVE KARMAN CANNULA MANUAL VACUUM ASPIRATOR**

Ref: 209

209 01 Single Valve Karman Cannula Manual Vacuum Aspirator Set (Aspirator + Cannulas)  
 209 04 Single Valve Karman Cannula Manual Vacuum Aspirator (No Cannulas)

- Includes extra O-ring, stopper ring and 2cc silicone oil.
- Stable vacuum: Each product's vacuum is individually tested.
- Portable.
- 60cc barrel capable of producing a vacuum of 609.6-660.4 mm/Hg.
- It is safe and effective.
- It is affordable.
- Appropriate for many different clinical settings.
- Provides high levels of patient and supplier satisfaction.
- Outpatient Procedure.
- It is single use device.
- Our Single Valve Karman Cannula Manual Vacuum Aspirator Set includes sterile Karman Cannulas sizes 4, 5 and 6.

**Alternative Names**

- Suction Aspirator
- Uterine Evacuation Service
- Aspirating Syringe
- Miscarriage Management Kit
- Gynecological Aspirator
- Aspirator Kit
- Manual Uterine Aspirator
- Endometrial Aspirator

**DOUBLE VALVE KARMAN CANNULA MANUAL VACUUM ASPIRATOR**



Ref: 209

209 02 Double Valve Karman Cannula Manual Vacuum Aspirator Set (Aspirator + Cannulas)  
 209 05 Double Valve Karman Cannula Manual Vacuum Aspirator (No Cannulas)

- Includes extra O-ring, stopper ring and 2cc silicone oil.
- Adaptors included. Stable vacuum: each product's vacuum is individually tested.
- Portable.
- 60cc barrel capable of producing a vacuum of 609.6-660.4 mm/Hg.
- It is safe and effective.
- It is affordable.
- Appropriate for many different clinical settings.
- Provides high levels of patient and supplier satisfaction.
- Outpatient Procedure.
- It is single use device.
- Our Double Valve Karman Cannula Manual Vacuum Aspirator Set includes sterile Karman Cannulas sizes 6, 7, 8, 9, 10 and 12.

**Alternative Names**

- Suction Aspirator
- Uterine Evacuation Service
- Aspirating Syringe
- Miscarriage Management Kit
- Gynecological Aspirator
- Aspirator Kit
- Manual Uterine Aspirator
- Endometrial Aspirator



**KARMAN CANNULA**

Ref: 208

208 04 Karman Cannula No: 4    208 05 Karman Cannula No: 5    208 06 Karman Cannula No: 6    208 07 Karman Cannula No: 7  
 208 08 Karman Cannula No: 8    208 09 Karman Cannula No: 9    208 10 Karman Cannula No: 10    208 12 Karman Cannula No: 12

- Total Length: 22 cm.
- It has smooth cannula surface.
- It has rounded tip.
- It is packaged separately in a sterile pouch.
- It has the same color codes as the adapter, allowing easy matching.

**Alternative Names**

- MVA Cannula
- Intra-uterine Suction Cannula



**SURGICAL DRAPES**

**C-ARM SCOPY COVER SET (WITH CLIPS)**

Ref: 227

227 04 C-arm Scopy Cover Set with Clips (Large)    227 08 C-arm Scopy Cover Set with Clips (Medium)    227 11 C-arm Scopy Cover Set with Clips (Small)

- Includes two scopy covers and one C-arm cover with clips.
- The arm of the machine is easily covered by the cover with clips and prevents contamination.
- Stabilizer clips are included.
- It completely surrounds the equipment and prevents contamination.
- Compatible with Siemens, Philips, Ziehm and GE C-Arm machines.

**MINI C-ARM SCOPY COVER**

Ref: 227

227 33 Mini C-Arm Scopy Cover 137x160 cm (Small)    227 34 Mini C-Arm Scopy Cover 137x216 cm (Large)

- It completely surrounds the equipment and prevents contamination.
- Special U-shape design conforms easily to the device.
- Three durable belts are available to properly secure the device to the machine.
- Compatible with all mini c-arm models.





**C-ARM UNIVERSAL PROTECTION COVER STRAIGHT** Ref: 227/291

227 09 C-Arm Universal Protection Cover Straight 90x225 cm (Large)  
227 10 C-Arm Universal Protection Cover Straight 90x300 cm (XL)

- Specially folded for easy application.
- Directional labels available for correct application.
- Two durable belts are available to properly secure the device to the machine.
- More than 8 different sizes are available.
- It can be produced in different sizes upon request.

**Alternative Names** • C-arm Drape



**C-ARM UNIVERSAL PROTECTION COVER SLOPPED** Ref: 227

227 29 C-Arm Universal Protection Cover Slopped 75+42+200 cm  
227 38 C-Arm Universal Protection Cover Slopped 75+42+220 cm

- Specially folded for easy application.
- It has slopped tip designed to better adapt to the machine.
- It has a thickness of 40 microns.
- Closed one side.
- Durable multiple belts are available to properly secure the device to the machine.
- It can be produced in different sizes upon request.

**Alternative Names** • C-Arm Scopy Drape • Fitted C-Arm Drape • Guesstet C-Arm Drape

**C-ARM SCOPY COVER SET (WITHOUT CLIPS)** Ref: 227

227 30 C-Arm Scopy Cover Set Without Clips 30x195 cm (Scopy Cover 65 cm x Scopy Cover 75 cm)  
227 31 C-Arm Scopy Cover Set Without Clips 30x195 cm (Fluoroscopy Cover 75 x 75 cm)  
227 36 C-Arm Scopy Cover Set Without Clips 30x208 cm (Scopy Cover 65 cm x Scopy Cover 75 cm)  
227 37 C-Arm Scopy Cover Set Without Clips 32x208 cm (Fluoroscopy Cover 80x75 cm x Fluoroscopy Cover 65x65 cm)  
227 35 C-Arm Scopy Cover Set Without Clips 40x190 cm (Fluoroscopy Cover 75x75)

- Specially folded for easy application.
- Upper and lower stickers for easy identification.
- It consists of 1 unit C arm cover and 2 unit Scopy/Fluoroscopy cover.
- It is available in 5 different sizes.
- It can be produced in different sizes upon request.



**FLUOROSCOPY COVER** Ref: 222 01

222 04 Fluoroscopy Cover 75x75 cm 222 05 Fluoroscopy Cover 75x90 cm 291 23 Fluoroscopy Cover 90x75 cm

- Ref: 222 01: 90x60 cm standard size (more than 20 different sizes are available.)
- Compatible with all devices.
- It can be produced in different sizes upon request.

**Alternative Names** • Rectangular Banded Bag



**SCOPY COVER** Ref: 227/291

291 10 Scopy Cover Ø 50 cm 291 11 Scopy Cover Ø 60 cm

- Ref: 227 01: Ø65 cm standard diameter (more than 15 different sizes available).
- Transparent: Does not inhibit the functioning of the device during use.
- Compatible with all devices
- It can be produced in different sizes upon request.

**Alternative Names** • Circular Banded Bag • Banded Bag • Brightness Amplifier Cover  
• Image Intensifier Cover • Bonnet Bag

**CARDBOARD CAMERA COVER** Ref: 202

202 15 Cardboard Camera Cover 14x140 cm 202 01 Cardboard Camera Cover 14x250 cm 202 25 Cardboard Camera Cover 14x300 cm  
202 17 Cardboard Camera Cover 18x250 cm 202 27 Cardboard Camera Cover 18x300 cm

- It has centralized hanger hole.
  - Telescopically folded.
  - Single layer thickness is 40 micron.
  - Eight layer folding allows the cover to pass over the camera easier and faster.
  - Compatible with all endoscopic cameras.
- It can be produced in different sizes upon request.

**Alternative Names** • Disposable Endoscopy Camera Cover • Sterile Endoscopy Camera Cover





**CARDBOARD CAMERA COVER ELASTIC TIP (BLUE)**

Ref: 202 07

202 28 Cardboard Camera Cover Elastic Tip 14x140 cm  
 202 37 Cardboard Camera Cover Elastic Tip 14x300 cm  
 202 30 Cardboard Camera Cover Elastic Tip 18x300 cm  
 202 07 Cardboard Camera Cover Elastic Tip 14x250 cm  
 202 20 Cardboard Camera Cover Elastic Tip 18x250 cm

- Specially designed elastic tip provides easy application.
- Telescopically folded.
- Single layer thickness is 40 micron.
- Eight layer folding allows the cover to pass over the camera easier and faster.
- Compatible with all camera models.
- It can be produced in different sizes upon request.

**Alternative Names** • Disposable Endoscopy Camera Cover • Sterile Endoscopy Camera Cover

**CARDBOARD CAMERA COVER ELASTIC TIP (WHITE)**

Ref: 202 07

202 07WE Cardboard Camera Cover Elastomer Tip WxL 14x250 cm  
 202 37WE Cardboard Camera Cover Elastomer Tip WxL 14x300 cm  
 202 20WE Cardboard Camera Cover Elastomer Tip WxL 18x250 cm  
 202 30WE Cardboard Camera Cover Elastomer Tip WxL 18x300 cm

- Specially designed elastic tip provides easy application.
- Telescopically folded.
- Single layer thickness is 40 micron.
- Eight layer folding allows the cover to pass over the camera easier and faster.
- Compatible with all endoscopic camera models.
- It can be produced in different sizes upon request.

**Alternative Names** • Disposable Endoscopy Camera Cover • Sterile Endoscopy Camera Cover

**ACCORDION FOLDED CAMERA COVER**

Ref: 202

202 21 Accordion Folded Camera Cover 14x250 cm  
 202 22 Accordion Folded Camera Cover 14x300 cm

- It has centralized hanger hall.
- It is folded as an accordion.
- Single layer thickness is 40 micron.
- Compatible with all camera models.
- Cardboard insert allows the cover to pass over the camera easier and faster.
- It can be produced in different sizes upon request.

**Alternative Names** • Concertina Folded Camera Cover



**CIRCLED CAMERA COVER**

Ref: 202 02

202 02 Circled Camera Cover 14x250 cm

- Specially folded for fast application.
- It has a larger center hole than any other camera cover on the market.
- It has adhesive tapes for easy fixing. (Extra tape is provided with each product.)
- Open and triangular cut tips are available.

**Alternative Names** • Disposable Endoscopy Ring Camera Cover

**CIRCLED CAMERA COVER ELASTIC TIP**

Ref: 202 09

202 09 Circled Camera Cover Elastic Tip Blue 14x250 cm  
 202 31 Circled Camera Cover Elastic Tip Blue 15x250 cm  
 202 09WE Circled Camera Cover Elastic Tip White 14x250 cm

- Specially designed elastic tip provides easy application.
- Specially folded for fast application.
- It has a larger center hole than any other camera cover on the market.
- It has adhesive tape for easy fixing.
- Blue and white elastic tips are available.

**Alternative Names** • Disposable Endoscopy Ring Camera Cover with Elastic Tip

**CLOSED SYSTEM CAMERA COVER**

Ref: 202 41

202 41 Cardboard Camera Cover Closed System 14x250 cm  
 202 44 Circled Camera Cover Closed System 15x250 cm

- It has centralized hanger hall.
- Telescopically folded.
- Cardboard or circled models are available.
- Eight layer folding allows the cover to pass over the camera easier and faster.
- Allows endoscopes to be easily changed during the procedure without disturbing the sterile field.
- It has an easy open, clamping system with spring-loaded locking mechanism.
- Reduces procedure time and turnaround.
- It can be produced in different sizes upon request.

**Alternative Names** • Closed System Coupler Camera Cover




**MICROSCOPE DRAPE**

Ref: 204

- 204 01 Microscope Drape (48 mm Straight Lens)
- 204 02 Microscope Drape (65 mm Slopped Lens)
- 204 03 Microscope Drape (68 mm Straight Lens)
- 204 07 Microscope Drape (65 mm Straight Lens)

- It has specially folded for easy application.
- Super transparent lens system is available.
- Adhesive belts secure the cover safely in place.
- Accessories are easy to use when the cover is in place.
- Compatible with Carl Zeiss, Leica and Moller-Wedel operating microscope models.
- More than 50 different sizes are available.
- It can be produced in different sizes upon request.



**BONNET MICROSCOPE DRAPE**

Ref: 293

- 293 03 Bonnet Microscope Drape 65 mm Slopped Lens 40x34 cm
- 293 04 Bonnet Microscope Drape 65 mm Slopped Lens 76x60 cm

- It has specially folded for easy application.
- Super transparent lens system is available.
- It has elasticized opening for easy application and placement.
- It has 70 micron thickness.
- If there is no need to cover the entire microscope, this product is more affordable than other models.
- It can be produced in different sizes upon request.

**MICROSCOPE HANDLE COVER**

Ref: 293

- 293 01 Microscope Handle Cover (Oval Tip) 2x: 14x31 cm
- 293 02 Microscope Handle Cover (Triangular Tip) 2x: 15x36 cm

- Ideal for sterile protection of operating microscope handles.
- Oval and triangular tip options are available.
- There are 2 products in a sterile pouch.
- Each covers have 2 units of velcro tape for easy application.
- Sold separately, or together with our Microscope Drapes.
- It can be produced in different sizes upon request.

**PE PROBE COVERS / PE PROBE COVER SETS**

Ref: 309

- 309 01 PE/309 01 PES Probe Cover 8x61 cm
- 309 02 PE/309 02 PES Probe Cover 10x61 cm
- 309 03 PE/309 03 PES Probe Cover 10x122 cm
- 309 04 PE/309 04 PES Probe Cover 13x61 cm
- 309 05 PE/309 05 PES Probe Cover 13x122 cm
- 309 06 PE/309 06 PES Probe Cover 13x244 cm
- 309 07 PE/309 07 PES Probe Cover 15x61 cm
- 309 08 PE/309 08 PES Probe Cover 15x122 cm
- 309 09 PE/309 09 PES Probe Cover 15x230 cm
- 309 10 PE/309 10 PES Probe Cover 18x30 cm
- 309 11 PE/309 11 PES Probe Cover 18x60 cm
- 309 12 PE/309 12 PES Probe Cover 18x90 cm
- 309 13 PE/309 13 PES Probe Cover 18x120 cm
- 309 14 PE/309 14 PES Probe Cover 18x250 cm

- Specially folded for easy application.
- It is produce of polyethylene material.
- Single layer thickness is 40 micron.
- Covers can be packaged separately of the sets.
- Package contains sterile ultrason gel, 3 unit elastic band, adhesive tapes.
- It can be produced in different sizes upon request.



**Alternative Names**

- Ultrasound Probe Cover

**LIGHT HANDLE COVER**


Ref: 308

- 308 01 Light Handle Cover (Short)
- 308 02 Light Handle Cover (Long)

- Consists of Low Density Polyethylene (LDPE) film and two interlocking plastic rings.
- Short (135 mm) and long (165 mm) versions are available.
- Compatible with all sizes of different operating room light handles.
- The inner ring has specially designed tooth to grip the light handle and prevent it from falling.

**Alternative Names**

- Operating Lamb Handle Cover





STERILE EO CE 1984

**ENDO SEAL BAG**

Ref: 279

279 01 Endo Seal Bag (Medium)      279 02 Endo Seal Bag (Large)

- It is patent protected.
- It consists of inner nylon for storing sterile endoscope and reversible carrying bag
- Transparent area is available to easily identify the content.
- Clean/Contaminated markings to easily distinguish as clean (Green) or dirty (Red).
- Completely contains the soiled equipment reducing the risk of spreading infection.
- Closing the bag prevents contamination.

**Alternative Names**

- Endoscope Transport Bag
- Endoscope Processing Bag

**ENDOSCOPY MOUTHPIECE**

Ref: 207

207 09 Endoscopy Mouthpiece With Strap and Tooth Protection Adult  
 207 10 Endoscopy Mouthpiece With Strap Pediatric  
 207 11 Endoscopy Mouthpiece Without Frame and strap with tooth protection  
 207 12 Endoscopy Mouthpiece With Oxygen Entry without Strap  
 207 13 Endoscopy Mouthpiece Without Frame and Tooth Protection  
 207 14 Endoscopy Mouthpiece Econobite  
 207 15 Endoscopy Mouthpiece With Oxygen Entry And Strap  
 207 16 Endoscopy Mouthpiece Without Tooth Protection With Strap  
 207 17 Endoscopy Mouthpiece Without Strap With Tooth Protection And Frame



CE LATEX

- Smooth edges and smooth passing surface are available.
- Silicone tooth protection.
- It has strap to secure comfortably in place during the procedure.
- Compatible with mouth anatomy.
- Conform for different sizes of endoscopes. (pediatric and adult versions are available.)
- Endoscopy Mouthpiece with Oxygen Entry allows oxygen to be given through the oxygen entrance.

**Alternative Names**

- Endoscopy Mouth Guard
- Endoscopy Bite Block



STERILE EO CE 1984

**MUCOUS ASPIRATOR**

Ref: 229

229 01 Mucous Aspirator 40 ml      229 02 Mucous Aspirator 20 ml      229 03 Mucous Aspirator 100 ml      229 04 Mucous Aspirator 15 ml

- Compatible with all suction catheters between 4-20 ch.
- Sticker for patient information included.
- An additional cap is provided free of charge.
- Sterilized disposable device with transient use.
- Our Mucous Aspirator products consists of a 40 ml graduated chamber with a double-entry cap, a conical connector and hose (Hose Length 12 cm).
- Consists of a chamber, a double-entry cap, a connector with a cap, a conical connector and two connection hoses.

**Alternative Names**

- Mucous Sucker
- Mucous Extractor
- Mucous Extraction
- Container Mucous Trap
- Mucous Suction Pump

**MUCOUS ASPIRATOR WITH HOSE**

Ref: 305 01



STERILE EO CE 1984

- After the procedure the container can be easily sealed via the use of the hose.
- Compatible with all suction catheters between 4-20 ch.
- Up to 40 ml. of aspiration capacity is available.
- It has clearly defined accurate scale.
- A sticker is included in the set to use for patient information.
- An additional cap is provided free of charge.

**Alternative Names**

- Mucous Extractor with Tube
- Suction Trap
- Lukens Trap
- Mucous Trap with Tube Mucus Specimen Trap

**ARTHROSCOPY SET**

Ref: 203

203 01 Arthroscopy Set with Manual Pressure Pump      203 02 Arthroscopy Set



STERILE EO CE 1984

- Y-type irrigation set.
- Drilling ends without air inlet
- Airproof incision ends.
- Latches.
- Y-connector.
- Air-freeing cap.

**Alternative Names**

- Arthroscopy Fluid Management Device (Irrigation Console Accessory)
- Cytoscopy Set





STERILE EO

**DRY STERILE SURGICAL HAND BRUSH**

Ref: 205 03

- Ergonomic brush design (48x79x34 mm).
- Easy open package.
- It has soft bristles to avoid skin irritation.
- With or without nail pick options are available
- Smooth and soft sponge absorbs up to 90 ml liquid.
- For external use only.

**Alternative Names**

- Disposable Surgical Brush
- Sterile Surgical Brush
- Sterile Dry Surgical Brush

**DISINFECTANT SURGICAL HAND BRUSH**

Ref: 205

205 01 %4 Chlorhexidine Impregnated Surgical Hand Brush  
 205 02 %7.5 Providone-Iodine Impregnated Surgical Hand Brush

- It is easy to use Ergonomic brush design.
- Easy open package.
- It has soft bristles to avoid skin irritation.
- With or without nail pick options are available.
- Smooth and soft sponge absorbs up to 90 ml of liquid.
- Chlorhexidine and Povidone-Iodine impregnated options are available.
- For external use only.

**Alternative Names**

- Impregnated Scrub Brush
- Surgical Brush
- Impregnated Surgical Brush
- Surgical Sponge



BIOCIDAL



CE 1984

**VALVE URINE BAG**

Ref: 225

225 01 Valve Urine Bag  
 225 03 Valve Urine Bag with discharge

- 1500 ml Scale: the amount of urine can be easily monitored.
- It has advanced valve system (no spills, no leakage).
- Discharge system is available for easy emptying. (Optional)
- Anatomically designed for use by both males and females.
- The special handle system allows for the bag to be hung anywhere easily.

**SPIROMETER FILTERED MOUTPIECE**


Ref: 273

273 01 SFT 33 mm	273 02 SFT 30mm	273 03 SFT 33 mm	273 04 SFT 30 mm with Latch
273 05 SFT 26 mm	273 06 SFT 26 mm with Latch	273 07 SFT 33 mm with Latch (Large)	273 08 SFT 30 mm with Latch (Large)
273 09 SFT 33 mm (Large)	273 10 SFT 30 mm (Large)		

- Made out of electrostatic material.
- Bacterial Filter Efficiency Test (BFE) 99.99%.
- Viral Filtration Efficiency Test (VFE) 99.99%.
- Easy to open packaging.
- Dead space is less than 50 ml.
- Manufactured without the use of chemical adhesives.
- With or without Latch options are available.

**Alternative Names**

- Spirometry Mouthpiece



CE 1984

**VALVE EMESIS BAG**

Ref: 228

228 05 Valve Emesis Bag Transparent  
 228 06 Valve Emesis Bag White

- It has a clean and smooth surface.
- It consists of a cardboard mouthpiece and nylon bag parts.
- There is product information and a calibrated measurement scale on the bag part.
- It is available in 2 models, the bag part being transparent and white.
- There is a valve system at the mouth to prevent liquid or odor leakage.



CE

