medbar®

Medbar Tıbbi Malzemeler Tur. San. ve Tic. A.Ş.
1142 Sk. No:35 Sarnıç Gaziemir - İzmir / Türkiye

M: 444 11 76
T: +90 232 281 60 03
F: +90 232 281 66 48
info@medbar.com.tr

www.medbar.com.tr



Experience, Quality and Technology

With 45 + years experience, we are constantly working toward better solutions.

medbar®

info@medbar.com.tr
www.medbar.com.tr



As of today, Medbar has more than 200 employees and exports to more than 70 countries in the world.

We aim at major projects with our vision of

"Enhancing the quality of life by maximizing our contribution to the healthcare process".



STERILE EO CE 1984

CYLINDRICAL IV FLOW CONTROLLER — Ref: 133

133 01 Cylindrical IV Flow Controller
133 02 Cylindrical IV Flow Controller Long (62cm)
133 03 Cylindrical IV Flow Controller Rotary Luer Lock
133 04 Cylindrical IV Flow Controller without Y-Port
133 05 Cylindrical IV Flow Controller Needle Free

- It works on the priciple of gravity.
- 5 to 250ml/h scale range is available.
- Different ranges between 5 and 250ml/h can be provided upon request.
- It has acurate and constant flow rate.
- Universal male and female luer connections are available.
- Needle and needle-free y-port injection options are available.
- Standard set length is 47cm.
- Different hose lengths are available.

Alternative Names • IV Flow Controller • IV Precision Flow Device • Controflow



STERILE EO CE 1984

IV FLOW CONTROLLER — Ref: 233

233 01 IV Flow Controller with Y-Port
233 02 IV Flow Controller Long
233 03 IV Flow Controller Rotary Luer Lock
233 04 IV Flow Controller without Y-Port
233 05 IV Flow Controller Needle Free

- It works on the priciple of gravity.
- 2 to 350 ml/hour scale range is available.
- It has accurate and constant flow rate.
- Universal male and female luer connections are available.
- Needle and needle-free y-port injection options are available.
- Standard hose length is 47 cm. (Different hose lengths are available.)

Alternative Names • IV Flow Regulator • IV Precision Flow Device • Controflow



STERILE EO CE 1984

EXTENSION LINE (250 PSI) — Ref: 294

294 01 Extension Line 30cm (250 PSI)
294 02 Extension Line 50cm (250 PSI)
294 03 Extension Line 60cm (250 PSI)
294 04 Extension Line 75 cm (250 PSI)
294 05 Extension Line 90 cm (250 PSI)
294 06 Extension Line 100 cm (250 PSI)
294 07 Extension Line 120cm (250 PSI)
294 08 Extension Line 150 cm (250 PSI)

- 2.70 mm-4.10 mm standard diameter.
- Universal male and female luer connectors are available.
- There are varieties with M/F and M/M connectors in sizes 30cm, 50cm, 60cm, 75cm, 90cm, 100cm, 120cm, 150cm.
- It is resistant to pressure up to 250 PSI (17 bar).



STERILE EO CE 1984

PRESSURE RESISTANT EXTENSION LINE (800 PSI) — Ref: 295

295 01 Pressure Resistant Extension Line 30 cm (800 PSI)
295 02 Pressure Resistant Extension Line 50 cm (800 PSI)
295 03 Pressure Resistant Extension Line 60 cm (800 PSI)
295 04 Pressure Resistant Extension Line 75 cm (800 PSI)
295 05 Pressure Resistant Extension Line 90 cm (800 PSI)
295 06 Pressure Resistant Extension Line 100 cm (800 PSI)
295 07 Pressure Resistant Extension Line 120 cm (800 PSI)
295 08 Pressure Resistant Extension Line 150 cm (800 PSI)

- 1.50 mm-3.00 mm standard diameter.
- Appropriate for all types of infusion treatment.
- Universal male and female luer connectors are available.
- There are varieties with M/F and M/M connectors in sizes 30cm, 50cm, 60cm, 75cm, 90cm, 100cm, 120cm, 150cm.
- It is resistant to pressure up to 800 PSI (55 bar).



STERILE EO CE 1984

PRESSURE RESISTANT EXTENSION LINE (1200 PSI) — Ref: 297

297 01 Pressure Resistant Extension Line 100 cm (1200 PSI)
297 03 Pressure Resistant Extension Line 120 cm (1200 PSI)
297 05 Pressure Resistant Extension Line 150 cm (1200 PSI)
297 07 Pressure Resistant Extension Line 180 cm (1200 PSI)
297 09 Pressure Resistant Extension Line 200 cm (1200 PSI)
297 11 Pressure Resistant Extension Line 250 cm (1200 PSI)
297 13 Pressure Resistant Extension Line 300 cm (1200 PSI)

- 1.50 mm-4.10 mm standard diameter.
- Appropriate for all types of infusion treatment.
- Universal male and female luer connectors are available.
- There are varieties with M/F and M/M connectors in sizes 100cm, 120cm, 150cm, 180cm, 200cm, 250cm, 300cm.
- It is resistant to pressure up to 1200 PSI (83 bar).

Alternative Names • Connecting Line • Extension Set • Connection Set



STERILE EO CE 1984

SINGLE VALVE KARMAN CANNULA MANUAL VACUUM ASPIRATOR — Ref: 209

209 01 Single Valve Karman Cannula Manual Vacuum Aspirator Set (Aspirator + Cannulas)
209 04 Single Valve Karman Cannula Manual Vacuum Aspirator (No Cannulas)

- Includes extra O-ring, stopper ring and 2cc silicone oil.
- Stable vacuum: Each product's vacuum is individually tested.
- Portable.
- 60cc barrel capable of producing a vacuum of 610-660 mm/Hg.
- It is safe and effective.
- It is affordable.
- Appropriate for many different clinical settings.
- Provides high levels of patient and supplier satisfaction.
- Outpatient Procedure.
- It is single use device.
- Our Single Valve Karman Cannula Manual Vacuum Aspirator Set includes sterile Karman Cannulas sizes 4, 5 and 6.

Alternative Names
- Suction Aspirator Uterine
- Evacuation Service
- Aspirating Syringe
- Miscarriage Management Kit
- Gynecological Aspirator
- Aspirator Kit
- Manual Uterine Aspirator
- Endometrial Aspirator



STERILE EO CE 1984

DOUBLE VALVE KARMAN CANNULA MANUAL VACUUM ASPIRATOR — Ref: 209

209 02 Double Valve Karman Cannula Manual Vacuum Aspirator Set (Aspirator + Cannulas)
209 05 Double Valve Karman Cannula Manual Vacuum Aspirator (No Cannulas)

- Includes extra O-ring, stopper ring and 2cc silicone oil.
- Stable vacuum: Each product's vacuum is individually tested.
- Portable.
- 60cc barrel capable of producing a vacuum of 610-660 mm/Hg.
- It is safe and effective.
- It is affordable.
- Appropriate for many different clinical settings.
- Provides high levels of patient and supplier satisfaction.
- Outpatient Procedure.
- It is single use device.
- Our Double Valve Karman Cannula Manual Vacuum Aspirator Set includes Adaptors and sterile Karman Cannulas sizes 6, 7, 8, 9, 10 and 12.

Alternative Names
- Suction Aspirator Uterine
- Evacuation Service
- Aspirating Syringe
- Miscarriage Management Kit
- Gynecological Aspirator
- Aspirator Kit
- Manual Uterine Aspirator
- Endometrial Aspirator



STERILE EO CE 1984

KARMAN CANNULA — Ref: 208

208 04 Karman Cannula No: 4 208 05 Karman Cannula No: 5 208 06 Karman Cannula No: 6 208 07 Karman Cannula No: 7
208 08 Karman Cannula No: 8 208 09 Karman Cannula No: 9 208 10 Karman Cannula No: 10 208 12 Karman Cannula No: 12

- Total Length: 22 cm.
- It has smooth cannula surface.
- It has rounded tip.
- It is packaged separately in a sterile pouch.
- It has the same color codes as the adapter, allowing easy matching.

Alternative Names
- MVA Cannula
- Intra-uterine Suction Cannula

SURGICAL DRAPES



STERILE EO CE 1984



STERILE EO CE 1984

C-ARM SCOPY COVER SET (WITH CLIPS) — Ref: 227

227 04 C-arm Scopy Cover Set with Clips (Large)
227 08 C-arm Scopy Cover Set with Clips (Medium)
227 11 C-arm Scopy Cover Set with Clips (Small)


- Includes two scopy covers and one C-arm cover with clips.
- The arm of the machine is easily covered by the cover with clips and prevents contamination.
- It completely surrounds the equipment and prevents contamination.
- Compatible with Siemens, Philips, Ziehm and GE C-Arm machines.

MINI C-ARM SCOPY COVER — Ref: 227

227 33 Mini C-Arm Scopy Cover 137x160 cm (Small)
227 34 Mini C-Arm Scopy Cover 137x216 cm (Large)

- It completely surrounds the equipment and prevents contamination.
- Special U-shape design conforms easily to the device.
- Three durable belts are available to properly secure the device to the
- machine. Compatible with all mini C-arm models.

Alternative Names
- C-Arm Drape
- Sterile C-Arm Drape
- Disposable C-Arm Cover



STERILE EO CE 1984

C-ARM UNIVERSAL PROTECTION COVER STRAIGHT — Ref: 227

227 09 C-Arm Universal Protection Cover Straight 90x225 cm (Large)
227 10 C-Arm Universal Protection Cover Straight 90x300 cm (XL)

- Specially folded for easy application.
- Directional labels available for correct application.
- Two durable belts are available to properly secure the device to the machine.
- More than 8 different sizes are available.
- It can be produced in different sizes upon request.

Alternative Names • C-Arm Drape • Universal C-Arm Cover



STERILE EO CE 1984

C-ARM SCOPY COVER SET (WITHOUT CLIPS) — Ref: 227

227 31 C Arm Scopy Cover Set Without Clips 32x195cm (75x75x2)
227 56 C Arm Scopy Cover Set Without Clips 32*207 (75*75, 75*75) cm
227 57 C Arm Scopy Cover Set Without Clips 32*207 (75*75, 55*75) cm
227 58 C Arm Scopy Cover Set Without Clips 32*207 (80*100, 75*75) cm
227 60 C Arm Scopy Cover Set Without Clips 32*207 (Diameter 65) cm

- Specially folded for easy application.
- Upper and lower stickers for easy identification.
- It consists of 1 unit C arm cover and 2 unit Scopy/Fluoroscopy cover.
- It is available in 5 different sizes.
- It can be produced in different sizes upon request.

Alternative Names • 3-Piece C-Arm Drape Set • C-Arm Fluoroscopy Drape Set



STERILE EO CE 1984

FLUOROSCOPY COVER — Ref: 222

222 01 Fluoroscopy Cover 90*60 cm
222 04 Fluoroscopy Cover 75*75 cm
222 05 Fluoroscopy Cover 75*90 cm
222 22 Fluoroscopy Cover 90*75 cm

- Ref: 222 01: 90x60 cm standard size. (more than 20 different sizes are available.)
- Compatible with all devices.
- It can be produced in different sizes upon request.

Alternative Names • Rectangular Banded Bag



STERILE EO CE 1984



STERILE EO CE 1984



STERILE EO CE 1984

SCOPY COVER — Ref: 291

291 03 Scopy Cover Ø 65 cm
291 10 Scopy Cover Ø 50 cm
291 11 Scopy Cover Ø 60 cm
291 12 Scopy Cover Ø 70 cm
291 13 Scopy Cover Ø 80 cm
291 14 Scopy Cover Ø 55 cm

- Ref. 291-03: Ø65 cm standard diameter (more than 20 different sizes available).
- Transparent: Does not inhibit the functioning of the device during use.
- Compatible with all devices.
- It can be produced in different sizes upon request.

Alternative Names • Circular Banded Bag • Banded Bag • Brightness Amplifier Cover • Image Intensifier Cover • Bonnet Bag

CAMERA COVER WITHOUT CARDBOARD — Ref: 202

202 06 Camera Cover Without Cardboard 14*250 cm
202 13 Camera Cover Without Cardboard 18*250 cm
202 35 Camera Cover Without Cardboard Blue Elastic Tip 18*250 cm

- Telescopically folded.
- Single layer thickness is 40 micron.
- Eight layer folding allows the cover to pass over the camera easier and faster.
- Compatible with all endoscopic cameras.
- It can be produced in different sizes upon request.

Alternative Names • Disposable Endoscopy Camera Cover • Sterile Endoscopy Camera Cover

CARDBOARD CAMERA COVER — Ref: 202

202 01 Cardboard Camera Cover 14*250 cm
202 10 Cardboard Camera Cover 15*250 cm
202 17 Cardboard Camera Cover 18*250 cm
202 25 Cardboard Camera Cover 14*300 cm
202 27 Cardboard Camera Cover 18*300 cm

- It has centralized hanger hole.
- Telescopically folded.
- Single layer thickness is 40 micron.
- Eight layer folding allows the cover to pass over the camera easier and faster.
- Compatible with all endoscopic cameras.
- It can be produced in different sizes upon request.

Alternative Names • Disposable Endoscopy Camera Cover • Sterile Endoscopy Camera Cover



STERILE EO CE 1984

CARDBOARD CAMERA COVER ELASTIC TIP (BLUE) — Ref: 202

202 07 Cardboard Camera Cover Blue Elastic Tip 14*250cm
202 19 Cardboard Camera Cover Blue Elastic Tip 15*250 cm
202 20 Cardboard Camera Cover Blue Elastic Tip 18*250 cm
202 30 Cardboard Camera Cover Blue Elastic Tip 18*300 cm

- Specially designed elastic tip provides easy application.
- Telescopically folded.
- It has centralized hanger hole.
- Single layer thickness is 40 micron.
- Eight layer folding allows the cover to pass over the camera easier and faster.
- Compatible with all camera models.
- It can be produced in different sizes upon request.

Alternative Names • Disposable Endoscopy Camera Cover • Sterile Endoscopy Camera Cover



STERILE EO CE 1984

ACCORDION FOLDED CAMERA COVER — Ref: 202

202 21 Accordion Folded Camera Cover 14*250 cm

- It has centralized hanger hall.
- It is folded as an accordion.
- Single layer thickness is 40 micron.
- Compatible with all camera models.
- Cardboard insert allows the cover to pass over the camera easier and faster.
- It can be produced in different sizes upon request.

Alternative Names • Concertina Folded Camera Cover



STERILE EO CE 1984

CIRCLED CAMERA COVER — Ref: 202

202 02 Circled Camera Cover 14*250 cm
202 51 Circled Camera Cover Perforetad Tip 14*250 cm

- Specially folded for fast application.
- It has a larger center hole than any other camera cover on the market.
- It has adhesive tapes for easy fixing. (Extra tape is provided with each product.)
- Open and triangular cut tips are available.

Alternative Names • Disposable Endoscopy Ring Camera Cover



STERILE EO CE 1984



STERILE EO CE 1984



STERILE EO CE 1984

CIRCLED CAMERA COVER ELASTIC TIP — Ref: 202

202 09 Circled Camera Cover Blue Elastic Tip 14*250 cm

- Specially designed elastic tip provides easy application.
- Specially folded for fast application.
- It has a larger center hole than any other camera cover on the market.
- It has adhesive tape for easy fixing.

Alternative Names • Disposable Endoscopy Ring Camera Cover with Elastic Tip

CLOSED SYSTEM CAMERA COVER — Ref: 202

202 41 Closed System Cardboard Camera Cover 14*250cm
202 43 Closed System Cardboard Camera Cover 18*250 cm
202 45 Closed System Circled Camera Cover 14*250cm

- It has centralized hanger hall.
- Telescopically folded.
- Cardboard or circled models are available.
- Eight layer folding allows the cover to pass over the camera easier and faster.
- Allows endoscopes to be easily changed during the procedure without disturbing the sterile field.
- It has an easy open, clamping system with spring-loaded locking mechanism.
- Reduces procedure time and turnaround.
- It can be produced in different sizes upon request.

Alternative Names • Closed System Endoscopic Camera Cover

LIGHT HANDLE COVER — Ref: 308

308 01 Light Handle Cover (Short)
308 02 Light Handle Cover (Long)

- Consists of Low Density Polyethylene (LDPE) film.
- It has two interlocking plastic rings.
- Short (135 mm) and long (165 mm) versions are available.
- Compatible with all sizes of different operating room light handles.
- Specially designed inner ring to facilitate gripping the surgical light handle.
- Grip teeth designed to help reduce slipping or falling off during use.

Alternative Names • Operating Lamb Handle Cover



STERILE EO CE 1984

MICROSCOPE DRAPE Ref: 204

204 01 Microscope Drape 48mm Straight Lens 122*225cm
204 06 Microscope Drape 65 mm Slopped Lens 117*305cmm
204 21 Microscope Drape 68mm Straight Lens 122*300cm
204 85 Microscope Drape 68mm Slopped Lens Elastomer 115*260 cm
204 86 Microscope Drape 65mm Slopped Lens Elastomer 122*305 cm
204 87 Microscope Drape 68mm Slopped Lens Elastomer 117*203 cm

- It has specially folded for easy application.
- Super transparent lens system is available.
- Adhesive belts secure the cover safely in place.
- Accessories are easy to use when the cover is in place.
- Compatible with Carl Zeiss, Leica and Moller-Wedel operating microscope models.
- More than 50 different sizes are available.
- It can be produced in different sizes upon request.

Alternative Names

- Surgical Microscope Cover
- Surgical Microscope Drape
- Microscope Protective Cover
- Sterile Surgical Microscope Drape



STERILE EO CE 1984

MICROSCOPE HANDLE COVER Ref: 293

293 06 Microscope Handle Cover Elasticated Oval Tip 15*36 cm

- Specially designed to cover surgical microscope handles.
- Two handle covers included in each sterile package.
- Two adhesive tapes provided for easy and secure application.
- Flexible use: Can be used independently or in combination with a surgical microscope drape.
- Customizable sizes available upon request.

Alternative Names

- Surgical Microscope Handle Cover
- Surgical Microscope Handle Protector
- Operating Microscope Handle Cover
- Microscope Handle Drape



STERILE EO CE 1984



CE NON STERILE

ENDOSCOPY BAG Ref: 279

279 01 Endoscopy Bag (Endoseal Bag) Medium
279 02 Endoscopy Bag (Endoseal Bag) Large

- It is patent protected.
- It consists of inner nylon for storing sterile endoscope and reversible carrying bag.
- Transparent area is available to easily identify the content.
- Clean/Contaminated markings to easily distinguish as clean (Green) or dirty (Red).
- Completely contains the soiled equipment reducing the risk of spreading infection.
- Closing the bag prevents contamination.

Alternative Names

- Endoscope Transport Bag
- Endoscope Processing Bag

ENDOSCOPY MOUTHPIECE (NON-STERILE) Ref: 207

207 09 Endoscopy Mouthpiece With Strap and Tooth Protection Framed Adult
207 10 Endoscopy Mouthpiece With Strap Pediatric
207 11 Endoscopy Mouthpiece Without Frame and Strap With Tooth Protection
207 12 Endoscopy Mouthpiece With Oxygen Entry without Strap
207 13 Endoscopy Mouthpiece Without Frame and Tooth Protection Without Strap
207 14 Endoscopy Mouthpiece Econobite
207 15 Endoscopy Mouthpiece With Oxygen Entry and Strap
207 16 Endoscopy Mouthpiece Without Tooth Protection With Strap
207 17 Endoscopy Mouthpiece Without Strap With Tooth Protection And Frame

- Smooth edges and transition surfaces for a comfortable fit.
- Ergonomic design suitable for intraoral use.
- Models available with silicone teeth protectors.
- It has strap to secure comfortably in place during the procedure.
- Compatible with endoscopes of different diameters.
- Pediatric and adult versions available.
- Optional oxygen entry available on selected models.
- Endoscopy Mouthpiece with Oxygen Entry allows oxygen to be given through the oxygen entrance.

Alternative Names

- Endoscopy Mouth Guard
- Endoscopy Bite Block



STERILE EO CE 1984

MUCUS COLLECTION CONTAINER — Ref: 229

229 01 Mucus Collection Container (40 ml-65 cm Hose)
229 02 Mucus Collection Container (25 ml-25 cm Hose)
229 03 Mucus Collection Container (100 ml-100 cm Hose)
229 04 Mucus Collection Container (15ml-25 cm with Single Falcon Tube)

- Available in options compatible with 4–20 Ch suction catheters.
- Dedicated labeling area for recording patient information.
- Graduated collection chamber for volume measurement.
- Additional cap included.
- Designed for short-term use; sterile and single-use.
- Available in 15 mL, 25 mL, 40 mL, and 100 mL capacities.
- Various tubing lengths available, ranging from 10 to 100 cm.
- Consists of a collection chamber, dual-port lid, capped connector, tapered connector, and two connecting tubes.

Alternative Names

- Mucus Aspirator
- Mucus Collector
- Mucus Extractor
- Mucus Suction Pump
- Mucus Suction Container



STERILE EO CE 1984


MUCUS COLLECTION CONTAINER WITH LARGE-BORE TUBING — Ref: 305

305 01 Mucus Collection Container (40 ml-13 cm Thick Hose)

- Compatible with 4–20 Ch suction catheters.
- 40 mL collection capacity.
- Graduated collection chamber for volume measurement.
- Dedicated labeling area for recording patient information.
- Designed to allow closure using the tubing after the procedure.
- Additional cap included.
- Large-bore tubing design.

Alternative Names

- Mucus Extractor with Tubing
- Mucus Collector with Tubing
- Mucus Collector
- Lukens Trap



STERILE EO CE 1984

ARTHROSCOPY SET — Ref: 203

203 01 Arthroscopy Set (Y-Tur Set) with Pump
203 02 Arthroscopy Set (Y-Tur Set)

- Y-Type Irrigation
- Non-Vented Piercing Spikes
- Clamps for Flow Control
- Equipped with a Y-connector.
- Air vent cap included.
- Available with a squeeze bulb option.

Alternative Names

- Arthroscopy Fluid Management Device (Irrigation Console Accessory)
- Cystoscopy Set



STERILE EO

STERILE DRY SURGICAL HAND BRUSH — Ref: 205

205 03 Sterile Dry Surgical Hand Brush With Nailpick

- Ergonomic brush design (48x79x34 mm).
- Easy open package.
- It has soft bristles to avoid skin irritation.
- With or without nail pick options are available
- Smooth and soft sponge absorbs up to 90 ml liquid.
- For external use only.

Alternative Names

- Disposable Surgical Hand Brush
- Sterile Dry Surgical Brush



BIOCIDAL

DISINFECTANT SURGICAL HAND BRUSH — Ref: 205

205 01 Surgical Hand Brush Impregnated with %4 Chlorhexidine Gluconate(Hibitanol) With Nailpick
205 02 Surgical Hand Brush Impregnated with %7,5 Povidone Iodine (Betadine) With Nailpick

- Ergonomic brush design.
- Easy open package.
- It has soft bristles to avoid skin irritation.
- With or without nail pick options are available.
- Smooth and soft sponge absorbs up to 90 ml of liquid.
- Chlorhexidine and Povidone-Iodine impregnated options are available.
- For external use only.

Alternative Names

- Impregnated Scrub Brush
- Impregnated Surgical Brush
- Surgical Brush
- Surgical Sponge



VALVE URINE BAG — Ref: 225

225 01 Valve Urine Bag
225 03 Valve Urine Bag (with discharge)

- 1500 mL collection capacity.
- Graduated measurement scale for visual monitoring of urine volume.
- Valve system designed to prevent leakage and backflow.
- Optional drainage tap for easy emptying.
- Suitable for both male and female patients.
- Integrated hanging handle for convenient positioning and suspension of the bag.

Alternative Names • Urine Collection Bag • Urine Drainage Bag



URINE METER — Ref: 288

288 01 Urine Meter

- 500 mL graduated urine measurement chamber.
- 2000 mL urine collection bag.
- Total system capacity of 2500 mL.
- Designed for monitoring hourly and periodic urine output.
- Easy-to-read graduated measurement scale.
- System designed to allow urine transfer from the measurement chamber
- to the collection bag. Detachable 2000 mL collection bag.
- Designed to help maintain a closed urinary drainage system.
- Hydrophobic Filter.
- Integrated hanging system for secure attachment to a patient bed or
- suitable support. Drainage tap for convenient urine emptying.
- Drainage tubing suitable for catheter connection.
- Sterile and single-use.

Alternative Names • Urinary Measurement System • Urine Output Monitoring System • Urometer System




CHEST DRAINAGE SET — Ref: 400

400 01 Chest Drainage Set (2000 ml)
400 02 Chest Drainage Set (700 ml)

- Available in 2000 mL and 700 mL collection capacities.
- Graduated fluid collection chamber.
- Easy-to-read measurement scale for visual monitoring of drainage volume.
- Compact and Easy-to-Use design.
- Designed for drainage of air and/or fluid from the pleural cavity.
- Designed to provide protection against positive pressure.
- Compatible with external vacuum systems.
- Integrated hanging system for secure attachment to a bed or stand.
- Sterile and single-use.

Alternative Names • Pleural Drainage System • Thoracic Drainage Unit • Closed Chest Drainage System



SPIROMETER FILTERED MOUTPIECE — Ref: 273

273 01 Spirometer Filtered Mouthpiece (33 mm) (Small)
273 02 Spirometer Filtered Mouthpiece (30 mm) (Small)
273 03 Spirometer Filtered Mouthpiece with Latch (33 mm) (Small)
273 04 Spirometer Filtered Mouthpiece with Latch (30 mm) (Small)
273 05 Spirometer Filtered Mouthpiece (26 mm) (Small)
273 06 Spirometer Filtered Mouthpiece with Latch (26 mm) (Small)

- Electrostatic filter design.
- 99.99% Bacterial Filtration Efficiency (BFE).
- 99.99% Viral Filtration Efficiency (VFE).
- Dead space of less than 50 mL.
- Manufactured without the use of chemical adhesives.
- Easy-open packaging design.
- Available with or without a latch.

Alternative Names • Spirometry Mouthpiece


medbar
Endustriyel Etilen Oksit Gaz Sterilizasyonu
TOUCH